

# XVIVO Organ Chamber™

Maintains the lung in a moist environment

The XVIVO Organ Chamber™ is a single-use, sterile disposable container intended to be used as a temporary receptacle for isolated donor lungs in preparation for eventual transplantation into a recipient.

## Features:

- Support the lungs and provides a reservoir to collect any fluids excreting from the organ
- Minimizes evaporation and maintains moisture within the chamber during EVLP
- Transparent lid allowing visual inspection
- Pre-cut openings for secure fastening of tubing
- Drain port and base designed to facilitate recycling of perfusate
- Opening for ventilation tubing
- Made of PETg and medical grade TPE. Not made with DEHP and latex.



For US: The XVIVO Organ Chamber is only for use with XPS™ and is included in the XPS Disposable Lung Kit™.

## Order information:

Reference no:

19020

Product:

1 x XVIVO Organ Chamber™

XVIVO Perfusion AB, Mlångans gata 10, SE-412 51 Göteborg, Sweden

Tel +46 31 788 21 50 | xvivogroup.com

XVIVO Organ Chamber™ is a trademark of XVIVO Perfusion AB.

© 2021 XVIVO Perfusion AB. All rights reserved.

Issued: 2021-10 Doc ID: MAR-16376-v.1.0

**XVIVO**